



Labor Dr. W. + U. Zollinger
Schärenmoosstrasse 105 8052 Zürich
Tel. 044 302 31 31 Fax 044 302 31 47
e-mail laborzollinger@laborzollinger.ch



STS
SCHWEIZERISCHER PRÜFSTELLENDIENST
SERVICE SUISSE D'ESSAI
SERVIZIO DI PROVA IN SVIZZERA
SWISS TESTING SERVICE
STS 148

Zurich, 21 April 2012

Bio-Strath AG
Postfach
8032 Zürich

Test report

The following analysis report corresponds solely to the analyzed sample(s). More precise indications in regard to the used test methods, respectively the uncertainty of measurement can be provided to the customer on request. The test report may only be photocopied as a whole.

Order No. 1688/1
Date of entry: 16 February 2012
Beginning of the analysis: 21 February 2012

Samples:
1688/1 Strath Aufbaupräparat 250ml, L 1111-1

Analytical method: Test of effectiveness of antimicrobial preservatives according to Ph. Eur. 7

Results

1688/1 Strath Aufbaupräparat 250ml, L 1111-1

-total count: 100
-moulds&yeasts: 10

Test organism	initial contamination cfu g ⁻¹	cfu g ⁻¹ in the product after				
		0 d	2 d	7 d	14 d	28 d
<i>Pseudomonas aeruginosa</i> ATCC 9027	1,1 Mio	< 10	< 10	< 10	< 10	
<i>Staphylococcus aureus</i> ATCC 6538	1,3 Mio	< 10	< 10	< 10	< 10	
<i>Escherichia coli</i> ATCC 8739	1,3 Mio	< 10	< 10	< 10	< 10	
<i>Candida albicans</i> ATCC 10231	220'000	350	< 10	< 10	< 10	
<i>Aspergillus brasiliensis</i> ATCC 16404	170'000	35'000	23'000	19'000	16'000	180

"cfu" means "colony forming units"
"< 10" means "less than 10"

Dr. W. Zollinger
CEO